

Protect-all

Innovative Laminating Solutions

Protecting You With

FACE SHIELD


Reorder Number: PA7001510SH

Help Stop the Spread with Lightweight Clear Face Shield


*Optically Clear
Elastic Band with Snap
Comfortable Foam Forehead Strip
Packed 12 per Bag*

Assembly Instructions:

Remove protective paper.
Attach elastic strap using snap on loose end
of the strap.
Place on head placing foam strip on
forehead.
Shield may be trimmed with scissors to
adjust length and shape.



*Protect-all makes no warranties, either express or implied,
that the face shield prevents infection or the transmission of
viruses or disease.*



Protect-all, Inc,

109 Badger Parkway
Darien, WI 53114
(262)724-3292
sales@protect-all.com